



Anaesthesia breathing systems

for GE Healthcare® Carestation®




Anaesthesia ▪ Breathing Systems

Anaesthesia breathing systems for GE Healthcare® Carestation®

A range of breathing systems designed specifically for anaesthesia Carestations from GE Healthcare®. The range includes four Flextube™ breathing systems, seven Silver Knight™ anti-microbial, three of which include a spirometry set, and one system manufactured from our expandable tubing, Compact II.

Breathing Systems for the ADU Carestation®



Code	Description	Box Qty.
2010085	22 mm Flextube breathing system for use with ADU Carestation® with fixed elbow, 2 L reservoir bag and bag limb, ≥1.6 m	15
2010109	22 mm Silver Knight anti-microbial Flextube breathing system for use with ADU Carestation® with fixed elbow, 2 L bag and limb, ≥1.6 m	12
2011008	22 mm Flextube breathing system for use with ADU Carestation® with fixed elbow, 2 L reservoir bag and bag limb, ≥2.4 m	12
2011112	22 mm Silver Knight anti-microbial Flextube breathing system for use with ADU Carestation® with fixed elbow, 2 L reservoir bag, limb and spirometry set, ≥2.4 m	12

Make an enquiry

Breathing systems for the Aisys®, Avance®, and Aespire® Carestation®

Code	Description	Box Qty.
2010088	22 mm Flextube breathing system for use with Aisys®, Avance® and Aespire® Carestation® with fixed elbow, 2 L reservoir bag and bag limb, ≥1.6 m	15
2010111	22 mm Silver Knight anti-microbial Flextube breathing system for use with Aisys®, Avance® and Aespire® Carestation® with fixed elbow, 2 L reservoir bag and limb, ≥1.6 m	15
2011106	22 mm Silver Knight anti-microbial Flextube breathing system for use with Aisys®, Avance® and Aespire® Carestation® with fixed elbow, 2 L reservoir bag and limb, ≥2.4 m	12
2011110	22 mm Silver Knight anti-microbial Flextube breathing system for use with Aisys®, Avance® and Aespire® Carestation® with fixed elbow, 2 L reservoir bag, limb and spirometry set, ≥2.4 m	12
2011012	22 mm Flextube breathing system for use with Aisys®, Avance® and Aespire® Carestation® with fixed elbow, 2 L reservoir bag and bag limb, ≥2.4 m	12

Make an enquiry

Illustrated Silver Knight

Carestation is a registered trademark of Datex-Ohmeda, Inc.

Datex-Ohmeda, Inc., is a General Electric company. GE Healthcare is a registered trademark of General Electric Company.

IS5.22 • Issue 7 08.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

Please think before you print
 Save energy and paper.
 If you must print this information sheet please print it double sided.