



Compact™

Extendable anaesthetic breathing systems for both adult and paediatric applications



Anaesthesia ▪ Anaesthetic Breathing Systems



Compact™ extendable anaesthetic breathing systems providing a versatile solution

Optimum compression

Designed for economy, the unique Compact extendable anaesthetic breathing system can be compressed to almost one quarter of its extended length. Supplied in compressed form increasing the number of systems per box, there can be significant space and cost savings compared with other breathing systems.

- More Compact systems per box
- Lower storage requirements
- Lower freight volume costs
- Improved disposal

Lightweight system

Lighter than other breathing systems ensuring that the 'pull' on the patient's airway device is kept to a minimum.

Adjustable length


Packed and supplied in the shortest length for simple extension to the length required for clinical use.


Easy positioning

The nature of the Compact system tube allows for it to be set in any position required. This results in multiple positioning possibilities, ideal for head and neck surgery where the system can be positioned away from the patient, and does not require any additional patient connections.



Adult breathing systems



Code	Description	Box Qty.
2150000	22 mm Compact breathing system with luer elbow, ≥ 1.5 m	70
2151000 (S*) 	22 mm Compact breathing system with luer elbow, ≥ 2.0 m	70 (20*)
2158000	22 mm Compact breathing system with luer elbow, ≥ 3.0 m	40
1526000	22 mm Compact limb ≥ 2.0 m	100


[Make an enquiry](#)

 Illustrated  Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 2151000S (sterile box quantity is shown in brackets).




Adult breathing systems with 3 L reservoir bag




Code	Description	Box Qty.
2152000	22 mm Compact breathing system with 3 L bag, luer elbow and limb, ≥ 2.0 m	35
2153000	22 mm Compact breathing system with 3 L bag, luer elbow, filter and limb, ≥ 2.0 m	35

Adult breathing systems with 2 L reservoir bag



Code	Description	Box Qty.
2154000 (S*)	22 mm Compact breathing system with 2 L bag, luer elbow and limb, ≥ 2.0 m	35 (20*)
2156000	22 mm Compact breathing system with 2 L bag, luer elbow and limb, ≥ 3.0 m	25
2157000	22 mm Compact breathing system with 2 L bag, luer elbow and limb, ≥ 3.0 m	25
2150001	22 mm Compact bagging system with 2 L bag, ≥ 1.5 m	70

Adult breathing system with water traps



Code	Description	Box Qty.
2155000	22 mm Compact breathing system with ported y-piece, water traps and limb, ≥ 3.0 m	12


Illustrated Sterile option available


(S*) Add an S to the seven digit code number for the sterile version of this product eg. 2154000S (sterile box quantity is shown in brackets).



Compact™ paediatric circle breathing systems


Paediatric breathing systems




Code	Description	Box Qty.
2161000 	15 mm Compact breathing system with luer elbow, ≥ 1.5 m	85
2162000	15 mm Compact breathing system with luer elbow, ≥ 2.0 m	70
2163000	15 mm Compact breathing system with luer elbow, ≥ 3.0 m	40

[Make an enquiry](#)

Paediatric breathing systems with 1 L reservoir bag



Code	Description	Box Qty.
2164000 	15 mm Compact breathing system with 1 L bag luer elbow and limb, ≥ 2.0 m	35
2165000	15 mm Compact breathing system with 1 L bag luer elbow and limb, ≥ 3.0 m	25

[Make an enquiry](#) Illustrated

IS5.4 • Issue 9 08.25



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • India • Ireland • Italy • Japan • Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA