



## Smoothbore catheter mounts

Clear and crush-resistant tubing



**Airway Management** ▪ Patient Connections

## Smoothbore catheter mounts


A full range of Smoothbore catheter mounts for patient connection providing good flow characteristics, lower tubing compliance and clear, crush-resistant tubing. Available with a variety of patient elbows and port cap options.

### A guide to our ports and caps

For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.

#### Luer ports

For secure attachment of gas sampling lines with luer lock fittings




#### 7.6mm port and flip top cap

For insertion of all standard 7.6mm probes or for suctioning



#### Double flip top cap with seal


Bronchoscopy and suction ports in one patient connection




<b>Code</b>	5120000 (S*)	5180000 (S*)	5250000 (S*)
<b>Box Qty.</b>	50 (85*)	50 (85*)	50 (50*)
<b>Length</b>	120mm	180mm	250mm
<b>Compressible Volume</b>	30ml	37ml	51ml

**S** Sterile option available


(S\*) Add an S to the seven digit code number for the sterile version of this product eg. 5120000S (sterile box quantity is shown in brackets).




<b>Code</b>	5180104 (S*)
<b>Box Qty.</b>	50 (85*)
<b>Compressible Volume</b>	25ml




<b>Code</b>	5180197 (S*)
<b>Box Qty.</b>	50 (85*)
<b>Compressible Volume</b>	35ml



<b>Code</b>	5180114 (S*)
<b>Box Qty.</b>	50 (85*)
<b>Compressible Volume</b>	36ml



<b>Code</b>	5180193 (S*)
<b>Box Qty.</b>	50 (85*)
<b>Compressible Volume</b>	35ml



<b>Code</b>	5180199 (S*)
<b>Box Qty.</b>	50 (85*)
<b>Compressible Volume</b>	37ml



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
 T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

**Please think before you print**  
 Save energy and paper.  
 If you must print this information sheet please print it double sided.