




Solus™ Curve laryngeal mask airway

The laryngeal mask with a preformed airway tube



Airway Management • Airway Devices

Solus™ Curve laryngeal mask airways

Solus Curve is our laryngeal mask for those who prefer the insertion characteristics of a device with a curved airway tube.

As with all our laryngeal masks, Solus Curve is single use, latex free and supplied sterile.

Essential information

Clearly visible for quick reference



Integral inflation line

Reduces the risk of damage whilst in use

High quality valve

To ensure continual cuff integrity

The curve

Anatomical curve to aid insertion

Low friction material

Atraumatic insertion and removal

Classic cuff shape

Renowned for reliability



▶ Video available at www.intersurgical.com



Code	Description	Size	Weight	Box Qty.
8005050	Solus Curve, large adult, laryngeal mask airway	5	70+kg	20
8004050	Solus Curve, medium adult, laryngeal mask airway	4	50–70kg	20
8003050	Solus Curve, small adult, laryngeal mask airway	3	30–50kg	20

Make an enquiry

Ⓢ All supplied sterile

IS6.12 • Issue 1 02.17



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com

UK • Ireland • France • Germany • Spain • Portugal
 Italy • Benelux • Sweden • Denmark • Lithuania • Russia
 Czech Republic • Turkey • South Africa • China • Japan
 Taiwan • Philippines • USA • Canada • Colombia • Australia



Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485

Please think before you print
 Save energy and paper.
 If you must print this e-mail
 please print it double sided